

TOPPAN

TEINNOVATIONS™

DUAL-OVENABLE CPET TRAYS AND BOWLS CATALOG



**Natural Color and Black Trays
Made From 30% Post-Consumer
Recycled CPET Material Now
Available to Order**

Why Natural CPET?

- NIR Detectable
- Ability to be captured in a MRF
- Wider range of end markets when recycled





CPET Dual-Ovenable Trays and Bowls

If you'd like a custom depth, size, color or material, contact us at info@teinnovations.com.

Catalog Key: ● Available in Black ○ Available in Natural

6.67" X 5.04" COLLECTION



Item 7589-100
1 Compartment
6.67 x 5.04 x 1.12 in.
11.7 fl. oz.
832/case, 9984/skid
●○



Item 7589-150
1 Compartment
6.67 x 5.04 x 1.5 in.
18.5 fl. oz.
704/case, 8448/skid
●○



Item 7589-175
1 Compartment
6.67 x 5.04 x 1.75 in.
21.1 fl. oz.
2704/case, 10816/skid
●○



Item 7589-212
1 Compartment
6.67 x 5.04 x 2.12 in.
24.2 fl. oz.
680/case, 8160/skid
●○



Item 7589-240
1 Compartment
6.67 x 5.04 x 2.4 in.
26.9 fl. oz.
672/case, 8064/skid
●○



Item 7589-125
1 Compartment
6.67 x 5.04 x 1.25 in.
15.4 fl. oz.
672/case, 8064/skid
●○



Item 7589-150-2C
2 Compartment, 60/40
6.67 x 5.04 x 1.5 in.
16.5 fl. oz.
680/case, 8160/skid
●○



Item 7589-187-2C
2 Compartment, 60/40
6.67 x 5.04 x 1.87 in.
20.4 fl. oz.
736/case, 8832/skid
●○

6.25" X 1.75" COLLECTION



Item 7036
1 Compartment
Raised Bottom
6.25 x 1.75 in.
17.8 fl. oz.
2622/case, 10488/skid
●○



Item 8712
1 Compartment
Flat Bottom
6.25 x 1.75 in.
18.8 fl. oz.
2622/case, 10488/skid
●○

6.30" X 1.77" COLLECTION



Item 1972
1 Compartment
Flat Bottom
6.30 x 1.77 in.
20.9 fl. oz.
2346/case, 9384/skid
●○



Item 1993
1 Compartment
Raised Bottom
6.30 x 1.77 in.
19.2 fl. oz.
2346/case, 9384/skid
●○

ELITE COLLECTION



Item 1994
1 Compartment
Flat Bottom
5.88 x 1.82 in.
18.1 fl. oz.
2604/case, 10416/skid
●○



Item 3357
1 Compartment
Flat Bottom, Ribbed
4 x 1.49 in.
6.1 fl. oz.
6612/case, 26448/skid
●○



Item 600785
1 Compartment
Flat Bottom
4.19 x 4.19 x 1.27 in.
6.7 fl. oz.
864/case, 23328/skid
●○



We manufacture our trays from PET resin to resist cracking and breaking under frozen conditions and to provide functional use as a package for reheating in both microwave and conventional ovens. PET trays provide a reliable sealing surface for a variety of films. PET trays are suitable for modified atmosphere packaging (MAP) applications with proper lidding film.



CPET Dual-Ovenable Trays and Bowls

If you'd like a custom depth, size, color or material, contact us at info@teinnovations.com.

Catalog Key: ● Available in Black ○ Available in Natural

8.5" X 6.5" COLLECTION



Item 6003
2 Compartment
8.5 x 6.5 x 1.5 in.
25.9 fl. oz.
2000/case, 8000/skid
●○



Item 7590-125
1 Compartment
8.5 x 6.5 x 1.31 in.
29.3 fl. oz.
368/case, 4416/skid
●○



Item 7590-150
1 Compartment
8.5 x 6.5 x 1.5 in.
31.8 fl. oz.
368/case, 4416/skid
●○



Item 7590-187
1 Compartment
8.5 x 6.5 x 1.87 in.
40.1 fl. oz.
368/case, 4416/skid
●○



Item 7590-325
1 Compartment
8.5 x 6.5 x 3.25 in.
67.8 fl. oz.
312/case, 3744/skid
●○



Item 7590-212
1 Compartment
8.5 x 6.5 x 2.12 in.
45.9 fl. oz.
384/case, 4608/skid
●○



Item 7590-400
1 Compartment
8.5 x 6.5 x 4.04 in.
78.7 fl. oz.
272/case, 3264/skid
●○



Item 7590-150-2C
2 Compartment
60/40
8.5 x 6.5 x 1.5 in.
30.3 fl. oz.
384/case, 4608/skid
●○



Item 7590-187-2C
2 Compartment, 60/40
8.5 x 6.5 x 1.87 in.
36.9 fl. oz.
384/case, 4608/skid
●○



Item 7590-250-2C
2 Compartment
8.5 x 6.5 x 2.5 in.
45.7 fl. oz.
376/case, 4512/skid
●○



Item 7590-150-3C
3 Compartment
8.5 x 6.5 x 1.5 in.
28.5 fl. oz.
368/case, 4416/skid
●○



Item 7590-187-3C
3 Compartment
8.5 x 6.5 x 1.87 in.
35.5 fl. oz.
368/case, 4416/skid
●○



Item 7590-212-3C
3 Compartment
8.5 x 6.5 x 2.12 in.
39.4 fl. oz.
368/case, 4416/skid
●○

8.5" X 6.38" COLLECTION



Item 8050
3 Compartment
8.5 x 6.38 x 1.46 in.
25.7 fl. oz.
540/case, 4320/skid
●○



Item 8051
3 Compartment
8.5 x 6.38 x 1.94 in.
29.6 fl. oz.
540/case, 4320/skid
●○



Item 8087
1 Compartment
10.25 x 8.125 x 2.08 in.
63.8 fl. oz.
640/case, 2560/skid
●○



Item 8901
1 Compartment
10.25 x 8.125 x 1.73 in.
54.4 fl. oz.
720/case, 2880/skid
●○



Item 7339
2 Compartment, 50/50
10.25 x 8.125 x 2.08 in.
43.8 fl. oz.
560/case, 2240/skid
●○

10.25" X 8.125" COLLECTION



CPET Dual-Ovenable Trays and Bowls

If you'd like a custom depth, size, color or material, contact us at info@teinnovations.com.

Catalog Key: ● Available in Black ○ Available in Natural

8.5" X 6.44" COLLECTION



Item 7049
1 Compartment
8.5 x 6.44 x 1.57 in.
28.2 fl. oz.
1760/case, 7040/skid
●○



Item 7059
2 Compartment
50/50
8.5 x 6.44 x 1.5 in.
28.2 fl. oz.
1600/case, 6400/skid
●○



Item 7049-125
1 Compartment
8.5 x 6.44 x 1.26 in.
25.6 fl. oz.
1760/case, 7040/skid
●○



Item 8861 A 18HW
1 Compartment
7.36 x 5.39 x 1.4 in.
21.3 fl. oz.
2205/case, 8820/skid
●○



Item 8866A
2 Compartment, 50/50
7.36 x 5.39 x 1.66 in.
22 fl. oz.
2142/case, 8568/skid
●○

12.44" X 10.37" COLLECTION



Item 7591-164
1 Compartment
12.44 x 10.37 x 1.64 in.
80 fl. oz.
140/case, 1680/skid
●○



Item 7591-208
1 Compartment
12.44 x 10.37 x 2.08 in.
108.6 fl. oz.
140/case, 1680/skid
●○



Item 7591-250
1 Compartment
12.44 x 10.37 x 2.54 in.
131 fl. oz.
140/case, 1680/skid
●○

OTHER SIZES



Item 7043M
1 Compartment
6.73 x 5 x 1.38 in.
16.1 fl. oz.
3510/case, 14040/skid
●○



Item 7522
1 Compartment
7.36 x 5.39 x 1.89 in.
27.3 fl. oz.
2205/case, 8820/skid
●○



Item 8070
1 Compartment
7.68 x 5.81 x 1.19 in.
17.7 fl. oz.
2160/case, 8640/skid
●○

TOPPAN is a market-leading provider of bowl and tray solutions for the frozen and refrigerated food industries. We offer a broad range of options for all your food processing and food service needs.

- Dual-ovenable CPET is freezer, microwave and conventional oven safe from -30°F to 400°F.
- We are an industry leader in innovation and sustainability.
- Our CPET trays can contain up to 30% post-consumer recycled resin depending on our customers' sustainability requirements and resin availability.
- Industry's best low-temperature breakage performance.
- Industry's best product quality and part consistency.



CPET Dual-Ovenable Trays and Bowls

If you'd like a custom depth, size, color or material, contact us at info@teinnovations.com.

Catalog Key: ● Available in Black ○ Available in Natural

OTHER SIZES



Item 8034A
7.88 x 1 in.
20.1 fl. oz.
1904/case, 7616/skid
●○



Item 8903
1 Compartment
8.96 x 6.96 in.
48.4 fl. oz.
1280/case, 5120/skid
●○



Item 9707BA
1 Compartment
6.69 x 4.71 x 1.5 in.
15.8 fl. oz.
3744/case, 14976/skid
●○



Item 9741BN
1 Compartment
9.6 x 8 x 1.54 in.
45.8 fl. oz.
200/case, 6000/skid
●○



Item 7152
1 Compartment
Flat Bottom
6.25 (5.00 i.d.) x 1.52 in.
10.1 fl. oz.
3072/case, 12288/skid
●○



Item 8800
1 Compartment
Flat Bottom
5.22 x 1.82 in.
14 fl. oz.
2752/case, 11008/skid
●○



Item 7442
12 Compartment
16.63 x 12.63 x 1.85 in.
88.8 fl. oz.
150/case, 2400/skid
●○



Item 8062
1 Compartment
7.75 x 5.86 x 1.87 in.
31.3 fl. oz.
1615/case, 6460/skid
●○



Item 8539
2 Compartment,
7.68 x 5.81 x 1.2 in.
17.5 fl. oz.
2064/case, 8256/skid
●○



Item 8689A
7.88 x 1.61 in.
32.8 fl. oz.
1540/case, 6160/skid
●○



Item 8865A
2 Compartment, 50/50
7.36 x 5.39 x 1.4 in.
18.7 fl. oz.
2205/case, 8820/skid
●○

CPET trays are functional over a broad temperature range from -30° Fahrenheit to 400° Fahrenheit. Our PET trays meet the requirements of 21 CFR 177.1630 established by the U.S. Food and Drug Administration.



HEATING INSTRUCTIONS

Conventional Oven

- The tray must be placed on baking sheet and positioned in the center of the oven during heating. Do not position tray on top rack or near a heating element.
- This is to provide additional support and stability of the tray to ensure the safety of the consumer.
- Placing the tray at the center of the oven allows for even heat distribution.
- The tray must be removed from oven on baking sheet and allowed to sit for 5 minutes at room temperature prior to handling.

Microwave

- Follow the manufacturer's microwave directions for specific food types.
- Allow to cool for 2 minutes prior to removing from oven.

**Natural Color and Black Trays
Made From 30% Post-Consumer
Recycled CPET Material Now
Available to Order**

Why Natural CPET?

- NIR Detectable
- Ability to be captured in a MRF
- Wider range of end markets when recycled



TEINNOVATIONS™

phone 630 393 0012
email info@teinnovations.com
online www.teinnovations.com